# ***Ear & Forehead Infrared Thermometer***

Model 09-342

**Please read this instruction manual completely before operating this unit.**

**INSTRUCTION MANUAL • ENGLISH & ESPAÑOL**


**STOP!**  
**PLEASE ENSURE YOU HAVE ALL OF**  
**THE FOLLOWING COMPONENTS**  
**BEFORE USING YOUR EAR & FOREHEAD**  
**INFRARED THERMOMETER**

Thermometer

2 AAA Batteries

Instruction Manual  
Quick Start Guide

**IF YOU ARE MISSING ANY PARTS,  
INCLUDING INSERTS OR INSTRUCTION  
MANUALS, DO NOT RETURN TO PLACE OF  
PURCHASE. CONTACT CUSTOMER CARE AT  
866-326-1313.**

---

# INDEX

---

Introduction.....	4
Body Temperature and Fever.....	5
Product Features .....	6
Device & Display Symbols .....	7
Battery Installation.....	8
Tips for Measuring Body Temperature.....	9
Setting Language, Volume, Measurement Mode, Date/Time.....	10-11
Measuring Temperature by Ear .....	12-13
Measuring Temperature by Forehead .....	14-15
Measuring Object/Liquid Temperature .....	16
Recalling & Clearing Readings.....	17
Care & Safety Information .....	18
Cleaning & Storage .....	19
Troubleshooting.....	20
Device & Label Symbols.....	21
Electromagnetic Compatibility .....	22-23
FCC Statement .....	24
Product Specifications.....	25
Warranty .....	26
Instruccionnes en Espanol .....	27-52

---

**Toll-Free Customer Care Help Line:**

**1-866-326-1313**

**Monday – Friday 8:30 a.m. – 4:30 p.m. CST**

Distributed by:  
Veridian Healthcare  
1175 Lakeside Drive  
Gurnee, IL 60031  
[www.veridianhealthcare.com](http://www.veridianhealthcare.com)  
Made in China

#93-1442 01/21  
©2021 Veridian Healthcare, LLC

---

# INTRODUCTION

---

Thank you for purchasing this Thermometer. Please read this instruction manual completely before using your unit to ensure safe use, accurate results and long life of the thermometer.

This product is intended for household use; this product is not intended to diagnose or treat any health problem or disease and should not be used as a substitute for the consultation and advice of a physician or other medical professional. Thermometers should only be used under the supervision of an adult.

**EAR THERMOMETERS** (also known as tympanic or aural thermometers), are a great alternative to oral, rectal and underarm thermometers. The ear is an excellent site for measuring temperature, as the part of the brain that regulates body temperature shares a blood supply with the eardrum. This thermometer provides results several times faster than the fastest oral thermometer, making it ideal for infants and fidgety children. It can even be used while the patient is sleeping. Unlike readings from oral thermometers, factors such as smoking, drinking, or speaking do not affect the measurements taken by ear thermometers.

**FOREHEAD THERMOMETERS** take the temperature of people and objects with minimal contact with the patient, making it less invasive and threatening than other thermometers, as well as potentially more sanitary. It can even be used when the patient is sleeping. This thermometer also measures the temperature of objects and liquids, making it useful for determining if a baby's milk or bathwater is ready.

All people and objects emit infrared energy. The forehead thermometer measures that energy, taking multiple readings quickly and converting those data to a single peak temperature, which appears on the display.

With proper care, your new thermometer will provide your entire family with many years of safe, accurate and comfortable temperature readings.

---


## **BODY TEMPERATURE AND FEVER**

---

Body temperature can vary from 97°F(36.1°C) and 99°F(37.2°C) and still be considered normal. The medically accepted 'normal' body temperature is 98.6°F(37°C). Body temperature is commonly lower upon waking than at any point during the rest of an individual's active hours.

A fever is defined as a body temperature that is elevated above the normal for that person. It is important to determine what is normal for an individual before determining if a fever is present. Tracking an individual's temperature on a consistent basis, at the same body site and at the same time every day will help determine an individual's normal temperature.

# PRODUCT FEATURES

## DEVICE & DISPLAY SYMBOLS

---

	Measurement in progress
	Degrees Fahrenheit Degrees Celsius
	Forehead measurement mode
	Ear Measurement Mode
	Object Measurement Mode
	Memory Symbol/Record Count
	Low Battery
	Sound On (Note: no icon if sound is OFF)

**DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.**

---

# BATTERY INSTALLATION

---


This device requires 2 AAA batteries, included. A small Phillips head screwdriver (not included) is necessary to remove the battery cover.

When the battery power becomes low, the low battery symbol  will appear on the display. The thermometer can still be used during this time, but the batteries should be replaced as soon as possible.

If the batteries expire completely, "Lo" will be displayed along with the low battery symbol. Batteries must be replaced before the thermometer can be used.

To install or replace the batteries, follow these steps:

1. Use a Phillips head screwdriver to carefully remove the screw and battery cover.
2. Slide the battery cover off.
3. Carefully remove the old batteries.
4. Insert the new batteries, ensuring that the positive and negative ends are aligned correctly.
5. Replace the battery cover and secure with screw.
6. Dispose of the expired batteries in accordance with local environmental policies.

---

## TIPS FOR MEASURING BODY TEMPERATURE

---

- Wait at least one minute between consecutive readings, to ensure accuracy.
- The patient should remain still while the reading is being taken.
- If “Lo,” “Hi,” “Err,” or another error message appears, see the Troubleshooting section of this manual for instructions on how to resolve the problem.
- Infrared forehead temperature readings are equivalent to oral temperature readings.
- Readings taken while asleep should not be compared directly to readings taken while awake, as body temperature while asleep is typically lower.
- Do not take body temperature readings within 30 minutes of being outdoors, exercising or bathing.
- Level of activity, time of day, and other factors may affect body temperature. Keep records of your individual body temperature to use as a reference when judging a fever.

---

# SETTING LANGUAGE, VOLUME, MEASUREMENT MODE, DATE/TIME

---

This device has a number of custom settings.

- Language —English or Spanish
- Volume —6 levels; level 1 being 'off' (no audio announcements)
- Measuring Unit— Fahrenheit, Celsius
- Date, Time, Year—real time measurement history

When batteries are initially installed or replaced, the features must be set to your user preference as instructed below.

If you do not go through the setting operation, the device will still operate but will default to English language, mid level volume, Fahrenheit mode, and January 1 at midnight settings.

1. Press the **POWER** button to turn the thermometer on.
2. Press and hold the **SET** button for 3+ seconds to enter the 'setting' operation sequence; the unit will automatically move to Setting the Language.


## SETTING THE LANGUAGE

3. The current language will be announced; press **SET** to alternate between English and Spanish.

L01 - English

L02 - Spanish

4. Press **MEM** to confirm and move to next setting.


## SETTING THE VOLUME

5. The current volume will be demonstrated by announcing 'ENGLISH' (or SPANISH) and the volume level will flash.

V01 — OFF, no volume heard

V02 thru V06 — lowest to highest level


6. Press **SET** to scroll through the volume levels and **MEM** to confirm and move to the next setting.

# SETTING LANGUAGE, VOLUME, MEASUREMENT MODE, DATE/TIME

## SETTING THE MEASUREMENT TEMPERATURE °F/°C

- The current mode will be indicated by announcing 'FAHRENHEIT' (or CELSIUS) and the °F or °C will flash.
- Press **SET** to alternate and select the desired mode, then **MEM** to confirm and move to the next setting.

## SETTING THE DATE / TIME

- The current year will flash on the display. Press **SET** to advance to the desired year. Press **MEM** to confirm.

Repeat the same steps to set the Month, Date, Hour, and Minute:

**SET** to advance to desired month; **MEM** to confirm

**SET** to advance to desired date; **MEM** to confirm

**SET** to advance to desired hour; **MEM** to confirm

**SET** to advance to desired minute; **MEM** to confirm

**NOTE:** This device utilizes 24 hr Military Time.

Flashing minutes example

After final **MEM** to confirm the minutes, the unit will beep once, the display will show OFF and the unit will turn off. You may now proceed to taking the desired measurement.

### NOTE:

You may turn OFF the speaker volume at any time during use, ideal for nighttime or sleeping patient readings.

While in-use, press and release the SET button. The speaker icon (SHOW ICON HERE OR ABOVE BY NOTE) on the lower right of the display indicates if the sound is on or off. If it is shown, the sound is ON; if it is not shown, the sound is OFF. Press SET again to toggle between ON and OFF.


# MEASURING TEMPERATURE BY EAR

## NOTE:

Proper placement and position of the probe is crucial for effective and accurate measurement at the ear.

The ear canal should be 'straight' in order for the sensor probe to work most effectively. Gently pull the ear slightly up and away from the head with the opposite hand that is holding the thermometer

- For children under one year or who cannot sit on their own, lay the child flat and turn the head sideways, then gently pull the ear and place the probe.
- For older children and adults, stand at the shoulder or elbow of the seated patient and pull back the ear and place the probe.

To measure body temperature at the EAR:

1. Remove the clear cap from the tip of the thermometer, Fig 1.
2. Press and release the **EAR / POWER** button to turn the thermometer on.
3. The display will illuminate and the unit will briefly run a self test, Fig 2.


Fig. 1

Fig. 2

# MEASURING TEMPERATURE BY EAR

- The unit will beep twice, then show the most recent results in memory, Fig 3, and is now ready to take a new measurement.
- Gently insert the probe into the ear canal until a slight resistance is felt.
- Press the **EAR** button for a count of one-second and release; you will hear an audible beep.  
Do not remove the thermometer nor release the button until you hear the beep.
- The display will show the measurement results, Fig. 4. If the audio is enabled, the results will also be announced aloud.
- The LED fever alert signal light will illuminate, indicating where the results fall in the fever alert range.


Fig. 3


Fig. 4

## FEVER ALERT SIGNAL LIGHT:

**GREEN LIGHT**— If the results below 99.5°F (37.5°C) indicating no fever is present

**AMBER (YELLOW) LIGHT** —If the results are between 99.5°F - 100.2°F (37.5°C - 37.9°C) a slight fever is present.

**RED LIGHT** —If the results are 100.3 °F (38.0°C) and above, a fever is considered present.

- After 60-seconds of non-use, the thermometer will automatically beep and shut off.

# MEASURING TEMPERATURE BY FOREHEAD

To measure body temperature at the FOREHEAD:

1. The clear cap should be in place during forehead measurements, Fig 1. The hole in the cap allows the sensor to work effectively, while the rounded surface of the cap helps with the 'swipe' action of taking a measurement at the forehead.
2. Press and release the **EAR / POWER** button to turn the thermometer on.
3. The display will illuminate and the unit will briefly run a self test, Fig 2.
4. The unit will beep twice, then show the most recent results in memory, Fig 3, and is now ready to take a new measurement.
5. While looking at the patient, place the probe tip on one temple. Press and hold the **HEAD** button while slowly sliding the thermometer across the forehead to the opposite temple while keeping contact on the skin.


Fig. 1


Fig. 2


Fig. 3

## NOTE:

Proper placement and position of the probe is crucial for effective and accurate measurement.

The forehead should be clean and clear of excessive lotion or oil.

The temple aligns with the eyes on each side.

# MEASURING TEMPERATURE BY FOREHEAD

- Release the button, the unit will beep, and remove the thermometer from the patient's head. The display will show the measurement results, Fig. 4. If the audio is enabled, the results will also be announced aloud.
- The LED fever alert signal light will illuminate, indicating where the results fall in the fever alert range.


Fig. 4

## FEVER ALERT SIGNAL LIGHT:

**GREEN LIGHT**— If the results below 99.5°F (37.5°C) indicating no fever is present

**AMBER (YELLOW) LIGHT** —If the results are between 99.5°F - 100.2°F (37.5°C - 37.9°C) a slight fever is present.

**RED LIGHT** —If the results are 100.3 °F (38.0°C) and above, a fever is considered present.

- After 60-seconds of non-use, the thermometer will automatically beep and shut off.

# MEASURING OBJECT/LIQUID TEMPERATURE

This thermometer can measure objects and liquids from 32°F to 212°F (0°C - 100°C). To measure the temperature of an object or liquid, follow these steps:

1. Press the **EAR / POWER** button to turn the thermometer on. The unit will run a self-test and all symbols on the display will momentarily appear, Fig. 1.
2. The unit will beep twice, then show the most recent results in memory.
3. Press **MEM** and **SET** simultaneously for 3-seconds to enter object measurement mode. The object mode symbol will appear on the display.
4. Aim the infrared lens at the object being measured, holding the thermometer approximately 1.0" away, Fig. 2.
5. Press the **EAR/POWER** button for a count of one-second and release; you will hear an audible beep
6. The display will show the measurement results, Fig. 3. If the audio is enabled, the results will also be announced aloud. The fever light indicator does not activate during object measurements.
7. To revert to body measurement mode, press and hold the **MEM** and **SET** buttons simultaneously again for 3-seconds until the 'object' symbol disappears. The thermometer is now ready for body measurement at the Ear or Forehead.
8. After 60-seconds of non-use, the thermometer will automatically beep and shut off.


Fig. 1

Fig. 2


Fig. 3

# RECALLING & CLEARING READINGS

This thermometer automatically stores the last 30 readings taken. To view the readings stored in memory, follow these steps:

1. Press the **EAR / POWER** button to turn the thermometer on.
2. Press the **MEM** button to recall the most recent measurement, Fig. 1.
3. The display will show the measurement results, date/time and measurement type – ear, forehead or object. If the audio is enabled, the results will be announced aloud.
4. Press and release the **MEM** button to scroll through the readings from newest to oldest. The index number in the bottom left of the display indicates the sequential number of the stored reading on the screen, Fig. 2.
5. After 30 measurements have been stored in the memory, the oldest reading will be automatically deleted.


Fig. 1

Fig. 1

After review, press the **EAR/POWER** button to enter measurement mode, or stop-using and the thermometer will automatically shut off after 60-seconds of non-use.

---

## CARE & SAFETY INFORMATION

---

Improper care or use of your thermometer may result in injury, damage to the unit or inaccurate readings, which can lead to a misdiagnosed fever. Following these instructions will ensure the thermometer's accuracy and long life.

- Never allow the thermometer to be used without adult supervision.
- The infrared technology used in this thermometer is safe and is not harmful in any way. However, we recommend avoiding looking into the lens during use.
- Clean the lens before and after each use.
- Do not touch the lens with fingers.
- Keep the protective clear cap in place when not in use.
- Dropping the thermometer or causing other heavy impact may damage it.
- Do not store the thermometer in direct sunlight, dusty or humid environments, or extreme temperatures.
- Do not disassemble the thermometer or attempt to repair it.
- If a reading cannot be taken, and the troubleshooting steps do not fix the problem, the thermometer must be replaced.

---




## CLEANING & STORAGE

---

The lens should be cleaned before and after each use, to eliminate dirt and fingerprints and ensure the accuracy of your readings. Proper storage will prevent degradation and damage to the unit.

### CLEANING:

- An alcohol swab, or a soft cotton cloth moistened with 70% isopropyl alcohol, should be used to gently cleanse the lens. Do not use water to wash the thermometer.
- The body of the thermometer can be cleaned with a soft, dry cloth.
- Never use a paper towel to clean the probe tip, as it may scratch the lens.
- Do not immerse the thermometer or allow liquid to enter the casing.
- Do not use cleaning agents other than alcohol to clean the thermometer; abrasives, thinners or gasoline will damage the unit.

### STORAGE:

- Store the unit in an area away from sunlight, dust and humidity.
- Protect the thermometer from mechanical shock or heavy impact.
- Avoid extended exposure to extreme temperatures. The thermometer should be stored in an environment between  $-13^{\circ}\text{F}$  and  $131^{\circ}\text{F}$  ( $-25^{\circ}\text{C}$ — $55^{\circ}\text{C}$ ), and with relative humidity of 95% or lower.
- Remove the batteries during extended storage.

# TROUBLE SHOOTING

Performance may be degraded should one or more of the following occur:

- Operation outside stated temperature and humidity range.
- Storage outside state temperature and humidity range.
- Mechanical shock.
- Patient temperature is below ambient temperature.

The following chart indicates error messages that may appear on your display, why they appear, and how you may solve the problem.

Error Symbol	Problem or Cause	Solution
	Low battery.	Replace batteries before using thermometer again.
 (blank display)	Unit turned off or battery problem.	Press Power button to activate unit. Check battery type or polarity (position in the unit). Replace battery.
	Temperature measurement is below typical temperature range.	Reposition the thermometer and try again.
	Temperature measurement is above typical temperature range.	Reposition the thermometer and try again.
	Ambient air temperature is outside of the acceptable operating range (60° F - 104° F or 16° C - 40° C).	Use the thermometer in an environment within the permissible temperature range.

# DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out its intended function; in the case of this device the cuff is the Type BF Applied Part
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal
	Manufacturer—This symbol indicates the manufacture

# ELECTROMAGNETIC COMPATIBILITY

09-342 series is intended for use in the electromagnetic environment specified below. The customer or the user of 09-342 series must make sure that it is used in such an environment.

## Guidance and manufacturer's declaration - Electromagnetic emissions

Phenomenon	Professional healthcare facility environment a)	HOME HEALTHCARE ENVIRONMENT a)
Conducted and radiated RF EMISSIONS	a)	CISPR 11 Group 1 Class B
Harmonic distortion	Not Applicable	
Voltage fluctuations and flickering	Not Applicable	
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.</p> <p>b) The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard.</p>		

# ELECTROMAGNETIC COMPATIBILITY

## Guidance and manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8kV contact ± 2 kV, ±4kV ±, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m b) 80MHz - 2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	<p>COMPLIANT</p> <p>NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the 09-342 series can be requested from Veridian Healthcare using the contact information provided in this manual. However, it is advisable to keep the electromechanical aerosol equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimise possible interference.</p>	
RATED power frequency magnetic fields.	IEC 61000-4-8	30 A/m c) 50 Hz or 60 Hz	
<p>a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.</p> <p>b) Before modulation is applied.</p> <p>c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields.</p>			

---

# FCC STATEMENT

---

## IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference; and (2) this device must accept any interference received, including interference that may cause undesirable operation.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

**CAUTION:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**CAUTION:** The equipment complies with FCC RF exposure limits set forth for an uncontrolled environment.

**CAUTION:** The equipment must not be co-located or operating in conjunction with any other antenna or transmitter.

---

# SPECIFICATIONS

---

<b>Name</b>	Ear & Forehead Infrared Thermometer
<b>Model Number</b>	09-342
<b>Display System</b>	Liquid Crystal Display with 0.1° F (0.1° C) resolution
<b>Measuring Range</b>	Human body: 93.2°F – 109.4°F (34° C – 43° C) Object: 32.0°F – 212.0°F (0° C – 100° C)
<b>Accuracy</b>	Human body: 93°F – 109°F (34°C – 43°C) +/- 0.4°F (0.2°C) Object: ±2.7°F (±1.5°C) <86°F (30°C); ±5% ≥86°F (30°C)
<b>Memory</b>	Automatically stores up to 30 readings.
<b>Tone</b>	Two short beeps sounds when the unit is turned on and when measurement is complete; audible announcements of all results if volume is level 2 and up
<b>Operating Temperature</b>	60° F to 104° F (16.0° C to 40.0° C) with relative humidity up to 95% (non condensing)
<b>Storage Temperature</b>	-13.0° F to 131° F (-25.0° C to 55.0° C) with relative humidity up to 95% (non condensing)
<b>Shut-Off</b>	Powered-off automatically after every 60 seconds of non-use
<b>Battery</b>	2 – 1.5v alkaline AAA size
<b>Unit Dimensions</b>	5-1/4" (L) x 1-1/2" (W) x 1-3/4" (H) (134.3mm x 42.4mm x 47.6mm)
<b>Weight</b>	4.27 ounces (65 grams) without battery
<b>Accessories</b>	Instruction Manual, Quick Start Guide, two AAA size batteries
<b>ASTM</b>	E1965–98

Specifications are subject to change without notice.

---

## THREE-YEAR LIMITED WARRANTY

---

The warrantor guarantees that its Thermometer will be free from manufacturing defects for a period of three years from the date of purchase. This warranty does not cover batteries or other power sources that may be provided with or used with the Infrared Thermometer. This warranty is voided if the Thermometer is misused or abused in any manner.

If the Thermometer fails to operate during the time the original purchaser owns it, return it postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include your name, address, phone number and brief description of the specific problem with your unit. Take care to carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. The warrantor will repair or replace the defective unit at our option. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

This warranty constitutes the warrantor's only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. The warrantor will not be responsible for any indirect, incidental, special consequential or punitive damages or other loss, including, but not limited to damage to or loss of other property or liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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