PRODIGY®

NO CODING BLOOD GLUCOSE TEST STRIPS

For use with PRODIGY AutoCode®, Pocket® and Voice® Blood Glucose Monitoring Systems.

Patient Cautions:

 For in vitro diagnostic (external) use only.

 For single patient use only.

 Please review these instructions and the blood glucose meters' Owner's Manual before using PRODIGY No Coding Blood Glucose Test Strips.

Critically ill patients should not use a blood glucose meter to monitor their glucose levels and should consult with their healthcare professional.

INTENDED USE

PRODIGY No Coding Blood Glucose Test Strips are intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples. PRODIGY No Coding test strips are intended to be used by a single person and should not be shared.

The PRODIGY Blood Glucose Monitoring System is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes and is an at-home aid to monitor the effectiveness of diabetes control. The PRODIGY Blood Glucose Monitoring System should not be used for the diagnosis or screening of diabetes, or for neonatal use.

STORAGE AND HANDLING

To ensure that your test strips are working effectively, please use them according to the following recommendations:

Important: Do not use the test strips after the expiration date, as the results may be inaccurate.

✎ Write the expiration date on the test strip vial when it is first opened. Test strips expire 90 days after opening.

 Store the test strips in a dry place between 39.2°–104°F (4°–40°C) and between 10–85% relative humidity. Keep the test strips away from direct sunlight.

✎ Test strips must be stored in their original vial ONLY. Do not transfer them to a new vial or other container.

✎ Do not touch the test strips when your hands are wet.

✎ Close the test strip's vial lid immediately after removing a strip and use each strip promptly.

✎ KEEP THE VIAL LID CLOSED WHEN NOT IN USE.

✎ Do not bend, cut or alter the test strip. Doing so will lead to inaccurate results.

HOW IT WORKS

PRODIGY No Coding test strips provide a quantitative measurement of glucose in whole blood from the fingertip, palm, forearm, upper arm, calf or thigh. The blood is drawn into the absorbent channel of the test strip through capillary action, and the blood glucose level is measured by the meter. Results are displayed on the PRODIGY blood glucose meter. Only 0.7 µL of blood is required and the reaction time is 7 seconds.

MATERIALS PROVIDED

PRODIGY No Coding Blood Glucose Test Strips

MATERIALS REQUIRED BUT NOT PROVIDED

PRODIGY Blood Glucose Meter

PRODIGY Lancets

PRODIGY Control Solution

Please purchase these items separately.

FOR ACCURATE RESULTS

PRODIGY No Coding test strips should be used with a fresh blood sample. For best results, it is recommended that you read the following procedure and the glucose meter's Owner's Manual before you obtain a blood sample. To get a drop of blood, follow these steps:

1. Wash your hands and the puncture site thoroughly with soap and water. Rinse and dry them well. Do not clean the testing site with alcohol.

2. Lance the puncture site to obtain a drop of blood. You may massage or gently rub the area to get the required amount of blood. Do not squeeze the puncture site excessively.

Important: A blood sample can also be obtained from the palm, forearm, upper arm, calf or thigh. Alternate Site Testing (AST) should be done only during steady-state times (when glucose is not changing rapidly). Different testing sites (AST) may lead to different results. Consult with your healthcare professional before testing with a blood sample obtained from AST locations.

Use AST only:

1. In a pre-meal or fasting state (more than 2 hours after the last meal).

2. Two hours or more after taking insulin.

3. Two hours or more after exercise.

4. During steady state of blood glucose conditions.

Do NOT use AST if:

1. You have reason to believe you have hypoglycemia or hyperglycemia.

2. Your routine glucose results are often fluctuating.

3. The test result is to be used to calibrate continuous glucose monitors (CGMs).

4. The test result is to be used for insulin dose calculations.

NO CODING REQUIRED TEST STRIPS

No calibration needed. Simply insert the test strip and proceed to the "CONTROL SOLUTION TESTING" or "BLOOD GLUCOSE TESTING" section.

CONTROL SOLUTION TESTING

When to perform a control solution test:

• When you begin using a new vial of test strips.

• At least once a week.

• Whenever you suspect that the meter or test strips are not working properly.

• When your blood glucose test results are not consistent with how you feel, or when you think your test results are not accurate.

• If you drop the meter.

• When you change the batteries.

PRODIGY Control Solution is used to check the testing technique and system performance. PRODIGY Blood Glucose Monitoring Systems use both PRODIGY Control Solution High and Low. When the control solution test result fall within the control range given for the lot of test strips being used, the system is working as it should and your testing technique is good. The control range can be found on the test strip vial label. Read your PRODIGY Blood Glucose Monitoring System's Owner's Manual for complete testing details.

Source: American Diabetes Association (September 22, 2014, <http://www.diabetes.org/diabetes-basics/diagnosis/>)

Questionable or Inconsistent Results:

If you are receiving test results that are unusual or inconsistent with how you are feeling:

• Make sure that the drop of blood completely fills the absorbent channel of the test strip.

• Confirm that the test strips are not past their expiration date.

BLOOD GLUCOSE TESTING

1. Fully insert the test strip into the meter with the contact bar entering the test strip slot. When the strip will go no further, the meter will turn on and perform several self-checks.

2. Obtain a blood sample of at least 0.7 µL. Hold the drop of blood to the absorbent channel of the test strip until it completely fills the confirmation window. The meter will beep and begin to count down.

3. After 7 seconds, the blood glucose test result will be displayed. The reading is automatically stored in the meter's memory. Shut the meter off by removing the test strip. Properly discard the used test strip.

Test Results:

Your blood glucose test results are displayed either in milligrams per deciliter (mg/dL, U.S. standard) or millimoles per liter (mmol/L).

The meter is able to display test results from 20 to 600 mg/dL or from 1.1 to 33.3 mmol/L. If your test result is below the minimum value, repeat the test with a new strip to confirm this reading. A reading lower than the minimum value indicates very low blood glucose levels, or severe hypoglycemia. If a low result is confirmed, seek medical attention immediately.

If your test result is above the maximum value, repeat the test with a new strip to confirm this reading. A reading higher than the maximum value indicates very high blood glucose levels, or severe hyperglycemia. If a high result is confirmed, seek medical attention immediately.

Reference Values:
The expected fasting blood glucose range for an adult person WITHOUT diabetes is less than 100 mg/dL (5.55 mmol/L) and less than 140 mg/dL (7.77 mmol/L) up to 2 hours after meals. Please consult your healthcare professional to establish your specific target ranges.

Source: American Diabetes Association (September 22, 2014, <http://www.diabetes.org/diabetes-basics/diagnosis/>)

Important: The control range can change with each new vial of test strips. Always use the control solution range on the label of your current vial of PRODIGY No Coding test strips.

• Check the performance of the meter and test strips by performing a control solution test. High or low blood glucose levels can indicate a possible serious medical condition. If you continue to get results that are unusually high or low, consult your healthcare professional.

Chemical components in sensor:

1. *Glucose oxidase (Aspergillus niger)* 0.4%.

2. *Potassium ferricyanide* 3.8%.

3. Non-reactive ingredients 95.8%.

TEST PRINCIPLE

The blood glucose test is based on the measurement of an electrical current generated by the reaction of glucose with the reagent of the test strip. The meter measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Limitations:

• Hematocrit: Hematocrit levels below 20% or above 60% can cause false results. Please consult your healthcare professional if you do not know your hematocrit level.

• Neonatal Use: The test strips should not be used for the testing of newborns.

• Concentrations of Acetaminophen > 8.0 mg/dL, Ascorbic acid > 5.0 mg/dL, Aspirin > 60 mg/dL, Bilirubin > 90 mg/dL, Cholesterol > 500mg/dL, Creatinine > 5.0 mg/dL, Dopamine > 2.0 mg/dL, Galactose > 900 mg/dL, Genticic acid > 5.0 mg/dL, Glutathione > 53 mg/dL, Hydroxyurea > 3.0 mg/dL, L-dopa > 10 mg/dL, Maltose > 900 mg/dL, Methylodopa > 3.0 mg/dL, Tolazamide > 100 mg/dL, Tolbutamide > 400 mg/dL, Triglyceride > 2,000 mg/dL, Uric acid > 8.0 mg/dL, and Xylose > 100 mg/dL may cause inaccurate results.

• The system exhibits interference from Acetaminophen. Do not use during Acetaminophen treatment.

• Use only fresh capillary whole blood.

• Altitude Effects: Test strips may be used at altitudes up to 10,742 feet (3,275 m) without affecting test results.

PERFORMANCE CHARACTERISTICS

The following results apply to PRODIGY No Coding Blood Glucose Test Strips used with PRODIGY blood glucose meters and compared with laboratory results.

Accuracy: Within ± 15 mg/dL (0.83 mmol/L) at glucose concentration < 75 mg/dL (4.2 mmol/L) and within ± 20% at glucose concentration ≥ 75 mg/dL (4.2 mmol/L).

Table 1 represents samples for glucose results less than (<) 75 mg/dL.

Number (%)	Within ±5 mg/dL	Within ±10 mg/dL	Within ±15 mg/dL
Finger	9/14 (64%)	13/14 (93%)	14/14 (100%)
Palm	5/14 (36%)	14/14 (100%)	14/14 (100%)
Forearm	3/14 (21%)	13/14 (93%)	14/14 (100%)
Upper arm	8/14 (57%)	14/14 (100%)	14/14 (100%)
Calf	6/14 (43%)	12/14 (86%)	14/14 (100%)
Thigh	4/14 (29%)	13/14 (93%)	14/14 (100%)

Table 2 represents samples for glucose results greater than or equal to (≥) 75 mg/dL.

Number (%)	Within ±5 %	Within ±10 %	Within ±15 %	Within ±20 %
Finger	30/86 (35%)	59/86 (69%)	85/86 (99%)	86/86 (100%)
Palm	29/86 (34%)	54/86 (63%)	83/86 (97%)	86/86 (100%)
Forearm	28/86 (33%)	59/86 (69%)	83/86 (97%)	86/86 (100%)
Upper arm	32/86 (37%)	63/86 (73%)	82/86 (95%)	86/86 (100%)
Calf	30/86 (35%)	53/86 (62%)	84/86 (98%)	86/86 (100%)
Thigh	23/86 (27%)	48/86 (56%)	83/86 (97%)	86/86 (100%)

Linear regression analysis:

Finger: $y = 1.0177x - 2.6955$; $R^2 = 0.975$

Palm: $y = 1.0437x - 6.3933$; $R^2 = 0.9688$

Forearm: $y = 0.9628x + 5.9635$; $R^2 = 0.9675$

Upper arm: $y = 1.0175x - 3.9493$; $R^2 = 0.9681$

Calf: $y = 1.0058x - 2.4028$; $R^2 = 0.9684$

Thigh: $y = 1.0137x - 3.3484$; $R^2 = 0.9576$

Precision: CVs(%) of intermediate precision and repeatability was less than 5%.

Within Run Repeatability (n=100)

Lot I	Mean (mg/dL)	43.1	78.8	129.4	195.4	372.2
SD	2.2	2.6	3.7	4.8	8.9	
CV	–	3.3%	2.8%	2.4%	2.7%	
Lot II	Mean (mg/dL)	43.2	78.8	129.5	194.9	327.1
SD	2.4	2.5	3.4	5.0	8.8	
CV	–	3.2%	2.6%	2.6%	2.7%	

Lot III	Mean (mg/dL)	43.6	79.4	129.7	194.7	328.9
SD	2.2	2.6	3.5	5.4	8.8	
CV	–	3.3%	2.7%	2.6%	2.7%	

Between Run Repeatability

Mean (mg/dL)	50.1	118	250	468	118	250	50.1	119	251
SD	2.6	3.8	7.6	2.7	3.8	7.8	2.5	3.7	6.2
CV(%)	–	3.2	3.2	–	3.2	3.1	–	3.1	2.5

The device has been certified to meet the following standards:
98/79/EC, IEC 60601-1, IEC 61010-1, IEC 60601-1-2

Please refer to your PRODIGY Blood Glucose Monitoring System's Owner's Manual for more detailed directions and illustrations.

Manufactured for:
Prodigy Diabetes Care, LLC
PO Box 481928, Charlotte, NC 28269 USA

For Customer Care call: 1-800-243-2636
Monday - Friday 8:30am - 5:30pm (EST)
www.prodigymeter.com

©2015 Prodigy is a registered trademark of Prodigy Diabetes Care, LLC. All rights reserved.
18-024 PDC2015 TSI01 Rev. 1 08/18/2018

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©2015 Prodigy es una marca registrada de Prodigy Diabetes Care, LLC. Todos los derechos reservados.