

Exam Tables

McKesson exam tables feature a contemporary design and an easy to clean laminate surface. These comfortable tables are available in two wooden finishes—natural and dark cherry—and a variety of vinyl upholstery colors to match other furniture in your exam room. For added comfort and convenience, each table features a standard table paper holder and table paper cutter and an optional flat or standard adjustable back rest. The McKesson exam table construction use flat panel components with bolt-through, pre-threaded hardware for easy and secure assembly. The H-brace construction, scuff-resistant wood coating, and stain-resistant vinyl provide durability. Tables are available with or without a bottom shelf and in two different table widths—27 or 30 inches.

Features and Benefits

- Seamless rounded corners
- Easy to assemble
- Premium, stain resistant, woven, knit-backed vinyl comes in a variety of colors
- Wood topcoat designed to resist scuffs, stains and abrasions
- Standard paper holder and paper cutter included with certain models
- Flat or adjustable backrest available
- 400-pound capacity under normal use
- Firm 2-inch foam padding
- Available in 27- or 30-inch widths
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM


With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Base Color Choices



Natural

Dark Cherry

Upholstery Color Choices



Black

Desert Tan

Royal Blue

Gun Metal

Wedgwood

ORDERING INFORMATION

MFR #	DESCRIPTION	ADDITIONAL FEATURES	QTY
2060	Exam Table with Adjustable Backrest, Paper Holder and Cutter, 27"	Specify color by chart	1 ea/cs
2060F	Exam Table with Flat Back, Paper Holder and Cutter, 27"	Specify color by chart	1 ea/cs
2061	Exam Table with Adjustable Backrest, Paper Holder and Cutter, 30"	Specify color by chart	1 ea/cs
2061F	Exam Table with Flat Back, Paper Holder and Cutter, 30"	Specify color by chart	1 ea/cs
2063	Exam Table with Shelf, Adjustable Backrest, Paper Holder and Cutter, 27"	Specify color by chart	1 ea/cs
2063F	Exam Table with Shelf, Flat Back, Paper Holder and Cutter, 27"	Specify color by chart	1 ea/cs
2064	Exam Table with Shelf, Adjustable Backrest, Paper Holder and Cutter, 30"	Specify color by chart	1 ea/cs
2064F	Exam Table with Shelf, Flat Back, Paper Holder and Cutter, 30"	Specify color by chart	1 ea/cs